

**ATCC Medium: 1887 Trypticase Soy Yeast Extract Agar/Broth  
Medium**

**Agar Medium**

Trypticase Soy Broth.....30.0 g  
Yeast Extract.....3.0 g  
Agar.....15.0 g  
DI Water.....1000 ml

Final pH should be 7.0-7.2. Boil to dissolve agar. Autoclave at 121°C.

**Broth Medium**

Trypticase Soy Broth.....30.0 g  
Yeast Extract.....3.0 g  
DI Water.....1000 ml

Final pH should be 7.0-7.2. Autoclave at 121°C.