

**ATCC Medium: 1887 Trypticase Soy Yeast Extract Agar/Broth
Medium**

Agar Medium

Trypticase Soy Broth.....30.0 g
Yeast Extract.....3.0 g
Agar.....15.0 g
DI Water.....1000 ml

Final pH should be 7.0-7.2. Boil to dissolve agar. Autoclave at 121°C.

Broth Medium

Trypticase Soy Broth.....30.0 g
Yeast Extract.....3.0 g
DI Water.....1000 ml

Final pH should be 7.0-7.2. Autoclave at 121°C.