

**ATCC Medium: 2737 Leeming & Notman agar Modified MLNA media:**

Bacteriological peptone.....	10g
Glucose.....	10g
Yeast extract.....	2.0 g
Ox bile, desiccated.....	8.0g
Glycerol.....	10mL
Glycerol monostearate.....	0.5g
Tween 60.....	5.0mL
Olive oil.....	20mL
Agar.....	15g
Deionized water.....	1.0 L

Medium is sterilized by autoclaving at 110C for 20 minutes. For liquid media, omit the agar.