

1652

Beef Extract Broth

Beef extract	3.0	g
Peptone	5.0	g
NaCl	5.0	g
Agar*, if necessary	10.0	g
Distilled water	to 1000.0	ml

Adjust medium for final pH 7.1.

Autoclave at 121 C for 15 minutes.

Cool to approximately 50 C and aseptically add the following sterile solutions:

0.1M MgSO ₄	0.8	ml
0.1M MnSO ₄	2.0	ml
0.1M CaCl ₂	10.0	ml

*For top agar, add 6.0 g/L