

ATCC medium: 2034 Autoinducer bioassay medium (AB)

NaCl.....17.53 g
MgSO₄.....6.02 g
Casamino acids, vitamin-free.....2.0 g
Distilled water.....960.0 ml

Adjust to pH 7.0 with KOH. Autoclave at 121C for 15 minutes. Cool to room temperature and aseptically add:

1 M Potassium phosphate, pH 7.0...10.0 ml
0.1 M L-Arginine.....10.0 ml
50% Glycerol.....20.0 ml