

**ATCC medium: 1233 Revised anthranilic acid medium for ATCC 31622**

Na<sub>2</sub>HPO<sub>4</sub> .....6.0 g  
KH<sub>2</sub>PO<sub>4</sub> .....3.0 g  
NaCl.....0.5 g  
NH<sub>4</sub>Cl .....1.0 g  
Distilled water.....950.0 ml

Autoclave basal medium for 15 minutes at 121C. Cool and aseptically add the following sterile solutions:

0.01 M CaCl<sub>2</sub> .....10.0 ml  
0.1 M MgSO<sub>4</sub> .....10.0 ml  
20% Glucose.....25.0 ml  
1% Anthranilic acid (in ethanol)....5.0 ml