

**ATCC Medium: 1069 YPAD Medium for MAK mutants of *Saccharomyces***

**Complete Medium**

Adenine sulfate.....	0.4 g
Yeast Extract.....	10.0 g
Peptone.....	20.0 g
Agar.....	20.0 g
Glucose Solution (added aseptically).....	100 ml
DI Water.....	900 ml

Autoclave at 121°C. When medium has cooled add the following solution:

**Glucose Solution**

Glucose.....	20.0 g
DI Water.....	100 ml

Filter sterilize.