

ATCC medium: 838 *Spirillum* nitrogen-fixing medium

KH₂PO₄0.4 g
K₂HPO₄0.1 g
MgSO₄ · 7H₂O0.2 g
NaCl.....0.1 g
CaCl₂.....0.02 g
FeCl₃.....0.01 g
Na₂MoO₄ · 2H₂O.....0.002 g
Sodium malate.....5.0 g
Yeast extract.....0.05 g
Distilled water.....1.0 L

Adjust pH to 7.2-7.4.