

ATCC Medium: 2089 Nutrient Agar/Broth with Naladixic Acid

For Solid Media

Nutrient Agar.....23.0 g
NaCl.....5.0 g
Distilled water.....975.0 ml

Autoclave at 121 C for 15 minutes. Cool to 50 C and aseptically add the following filter-sterilized component:

Naladixic acid0.1 g
Distilled water..... 25.0 ml
Filter sterilize

For Broth

Nutrient Broth.....8.0 g
NaCl.....5.0 g
Distilled water.....975.0 ml

Autoclave at 121 C for 15 minutes. Cool to 50 C and aseptically add the following filter-sterilized component:

Naladixic acid0.1 g
Distilled water..... 25.0 ml
Filter sterilize