

**ATCC medium: 1074 GC Medium with defined supplements**

*Base:*

Agar.....5.0 g  
GC Agar Base (BD 212307)....36.0 g  
Distilled water.....1.0 L

Autoclave for 15 minutes at 121C. Cool to 50C in water bath.

*Defined Supplements:*

Glucose.....40.0 g  
Glutamine.....1.0 g  
Coccarboxylase.....2.0 mg  
Ferric nitrate.....0.05 g  
Distilled water.....100.0 ml

Dissolve supplements in distilled water. Filter-sterilize. Add 1% v/v to cooled (50C) GC medium base. GC Broth contains the same components as GC Agar without the agar and cornstarch.