

ATCC medium: 1963 R agar with catalase

Peptone (BD 211677).....10.0 g
Yeast Extract (BD 212750).....5.0 g
Malt Extract (BD 218630).....5.0 g
Casamino Acids (BD 223050).....5.0 g
Beef Extract Paste (BD 212610).....2.0 g
Glycerol.....2.0 g
Tween 80.....50.0 mg
MgSO₄ . 7H₂O1.0 g
Agar.....20.0 g
Catalase (filter-sterilized).....60.0 mg
Distilled water.....1.0 L

Adjust medium for final pH 7.2. Autoclave at 121C for 15 minutes. Cool to 50C or less and aseptically add sterile catalase.