

ATCC medium: 412 Methanol medium for *Achromobacter*

Methanol*	20.0 ml
NH ₄ Cl	5.0 g
KH ₂ PO ₄	2.0 g
NaCl	0.5 g
MgSO ₄	0.2 g
FeSO ₄	2.0 mg
MnCl ₂	2.0 mg
Yeast extract	0.2 g
Distilled water	1.0 L

*Filter-sterilized and added aseptically to sterile medium

Adjust pH to 7.0. Autoclave at 121C for 15 minutes and then aseptically add methanol.