

**ATCC medium: 627 *Aspergillus* medium**

Yeast extract.....1.0 g  
Peptone.....1.0 g  
NaNO<sub>3</sub>.....6.0 g  
Casamino Acids (BD 223050).....1.0 g  
Agar.....15.0 g  
Adenine.....0.15 g  
Vitamin Solution (see below).....10.0 ml  
Distilled water to.....1.0 L

Adjust pH to 6.0. Autoclave at 121C for 10 minutes.

*Vitamin Solution:*

Biotin.....0.01 g  
Pyridoxine hydrochloride.....0.01 g  
Thiamine . HCl.....0.01 g  
Riboflavin.....0.01 g  
p-Aminobenzoic acid.....0.01 g  
Nicotinic acid.....0.01 g  
Distilled water to.....100.0 ml

Autoclave at 121C for 10 minutes.