

**ATCC medium: 876 Casein medium**

NaCl.....250.0 g  
HY Case SF (casein hydrolysate) (Sigma C9386)..5.0 g  
Yeast extract.....5.0 g  
MgCl<sub>2</sub> . 6H<sub>2</sub>O .....20.0 g  
KCl.....2.0 g  
CaCl<sub>2</sub> . 2H<sub>2</sub>O .....0.2 g  
Agar.....20.0 g  
Distilled water to.....1.0 L

Adjust final pH to 7.6.  
Autoclave at 121C for 15 minutes.