

ATCC medium: 1294 Medium A

Yeast extract.....	10.0 g
D-Glucose.....	20.0 g
Biotin (filter-sterilized).....	1.0 mg
Calcium D-(+)-pantothenate (filter-sterilized).....	1.0 mg
Agar.....	20.0 g
Distilled water.....	1.0 L

Combine all ingredients except Biotin and Calcium D-(+)-pantothenate, adjust pH to 7.3, autoclave at 121C for 15 minutes, then aseptically add the filter-sterilized components.