



QUALITY CONTROL STRAINS FOR COMMERCIAL ASSAYS



ATCC® | Credible leads to Incredible™

Maintain your outstanding quality control programs

When it comes to product development, safety and quality are of the utmost importance. That's why as part of our commitment to raising the standards for credibility in the life sciences, ATCC provides the authenticated and fully characterized reference materials needed to maintain outstanding quality control programs.

ATCC cultures are meticulously characterized and preserved to maintain these materials in a manner that permits reproducibility of results across time and among laboratories around the world. With nearly a century's worth of expertise in biomaterials handling, storage, and authentication, you can trust that ATCC quality control strains are backed by polyphasic identification techniques, best practices in stock management, and robust quality systems.

Come to the source for credible cultures for your quality control needs.

Visit our website at www.atcc.org/qc to learn more.



TABLE 1. BD Diagnostic Systems

Sets and Organisms	ATCC® No.
BD CRYSTAL™ ANR QC SET	
<i>Bacteroides fragilis</i>	25285™
<i>Fusobacterium varium</i>	27725™
<i>Lactobacillus acidophilus</i>	314™
<i>Parabacteroides distasonis</i>	8503™
<i>Peptoniphilus asaccharolyticus</i>	29743™
BD CRYSTAL™ E/NF QC SET	
<i>Acinetobacter lwoffii</i>	17925™
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i>	35030™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	33495™ or 33495-MINI-PACK™
<i>Proteus vulgaris</i>	8427™
<i>Pseudomonas aeruginosa</i>	35032™
BD CRYSTAL™ GRAM-POSITIVE QC SET	
<i>Brevibacillus brevis</i>	8246™
<i>Enterococcus faecalis</i>	19433™
<i>Staphylococcus epidermidis</i>	12228™ or 12228-MINI-PACK™
<i>Staphylococcus xylosum</i>	35033™
<i>Streptococcus pyogenes</i>	19615™ or 19615-MINI-PACK™

TABLE 1. BD Diagnostic Systems

Sets and Organisms	ATCC® No.
BD CRYSTAL™ RAPID GRAM-POSITIVE QC SET	
<i>Brevibacillus brevis</i>	8246™
<i>Enterococcus faecalis</i>	19433™
<i>Staphylococcus epidermidis</i>	12228™ or 12228-MINI-PACK™
<i>Streptococcus pneumoniae</i>	6303™
<i>Streptococcus pyogenes</i>	19615™ or 19615-MINI-PACK™
BD CRYSTAL™ NH QC SET	
<i>Aggregatibacter aphrophilus</i>	19415™
<i>Haemophilus influenzae</i>	35056™
<i>Kingella denitrificans</i>	33394™
<i>Moraxella catarrhalis</i>	25240™
<i>Neisseria lactamica</i>	49142™
BD PHOENIX™ GRAM-NEGATIVE AST NMIC/ID QC SET	
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Escherichia coli</i>	35218™ or 35218-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	700603™ or 700603-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
BD PHOENIX™ GRAM-NEGATIVE ID QC SET	
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
BD PHOENIX™ GRAM-POSITIVE AST PMIC/ID QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Enterococcus faecalis</i>	51299™ or 51299-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	25923™ or 25923-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29213™ or 29213-MINI-PACK™
BD PHOENIX™ GRAM-POSITIVE ID QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29213™ or 29213-MINI-PACK™
BD PHOENIX™ STREPTOCOCCUS SMIC/ID QC SET	
<i>Streptococcus agalactiae</i>	13813™ or 13813-MINI-PACK™
<i>Streptococcus pneumoniae</i>	49619™ or 49619-MINI-PACK™
BD PHOENIX™ STREPTOCOCCUS SMIC QC SET	
<i>Streptococcus pneumoniae</i>	49619™ or 49619-MINI-PACK™
BD GENEOHM™ VANR ASSAY	
<i>Enterococcus faecalis</i> (Positive Control)	51299™ or 51299-MINI-PACK™
<i>Enterococcus faecalis</i> (Negative Control)	29212™ or 29212-MINI-PACK™
<i>Enterococcus gallinarum</i> (Negative Control)	700425™
BD GENEOHM™ MRSA ACP	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Positive Control)	43300™ or 43300-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Negative Control)	25923™ or 25923-MINI-PACK™
Alternative Strain: <i>Staphylococcus epidermidis</i>	14990™ or 14990-MINI-PACK™
BD GENEOHM™ STREPB	
<i>Streptococcus</i> sp. (Positive Control)	12973™
<i>Streptococcus bovis</i> (Negative Control)	33317™

TABLE 1. BD Diagnostic Systems

Sets and Organisms	ATCC® No.
BD GENEHM™ STAPHSR	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Positive MRSA)	43300™ or 43300-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Positive MSSA)	25923™ or 25923-MINI-PACK™
<i>Staphylococcus epidermidis</i> (Negative Control)	14990™ or 14990-MINI-PACK™
BD GENEHM™ CDIFF	
<i>Clostridioides difficile</i> (Positive Control)	43255™
<i>Clostridioides difficile</i> (Negative Control)	700057™
BD MAX™ GBS	
<i>Streptococcus agalactiae</i> (Positive Control)	BAA-22™ or BAA-22-MINI-PACK™
BD MAX™ MRSA	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Positive Control)	43300™ or 43300-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Negative Control)	25923™ or 25923-MINI-PACK™
BD MAX™ CDIFF	
<i>Clostridioides difficile</i> (Positive Control)	43255™
<i>Clostridioides difficile</i> (Negative Control)	700057™
BD MAX™ STAPHSR	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Positive Control)	43300™ or 43300-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Positive Control)	29213™ or 29213-MINI-PACK™
<i>Staphylococcus epidermidis</i> (Negative Control)	12228™ or 12228-MINI-PACK™
BD MAX™ MRSA XT	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (Positive Control)	43300™ or 43300-MINI-PACK™
<i>Staphylococcus epidermidis</i> (Negative Control)	12228™ or 12228-MINI-PACK™
BD MAX™ ENTERIC BACTERIAL PANEL	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium (Positive Control)	14028™ or 14028-MINI-PACK™
<i>Shigella sonnei</i> (Positive Control)	9290™
<i>Escherichia coli</i> , stx 1a (Positive Control)	43890™
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> (Positive Control)	33291™
BD MAX™ ENTERIC PARASITE PANEL	
<i>Giardia intestinalis</i> (Positive Control)	30888D™
<i>Cryptosporidium parvum</i> (Positive Control)	PRA-67DQ™
<i>Entamoeba histolytica</i> (Positive Control)	30459DQ™
BD MAX™ CT/GC/TV	
<i>Chlamydia trachomatis</i> serovar D (Positive Control)	VR-885™ or serovar H VR-879™
<i>Neisseria gonorrhoeae</i> (Positive Control)	19424™ or 19424-MINI-PACK™
<i>Trichomonas vaginalis</i> (Positive Control)	30001™
BD ENTEROPLURI-TEST KIT *NOT FOR HUMAN DIAGNOSTIC USE	
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	13883™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
<i>Proteus mirabilis</i>	25933™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium	14028™ or 14028-MINI-PACK™

TABLE 2. Beckman Coulter

Sets and Organisms	ATCC® No.
MICROSCAN® ESBL PLUS QC SET	
<i>Escherichia coli</i> (MIC)	25922™ or 25922-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> (ESBL)	700603™ or 700603-MINI-PACK™
<i>Pseudomonas aeruginosa</i> (MIC)	27853™ or 27853-MINI-PACK™
MICROSCAN® HNID QC SET	
<i>Neisseria lactamica</i> (ID)	49142™
<i>Haemophilus influenzae</i> (ID)	49144™
<i>Aggregatibacter aphrophilus</i> (ID)	49146™
<i>Moraxella catarrhalis</i> (ID)	49143™
<i>Gardnerella vaginalis</i> (ID)	49145™
<i>Neisseria mucosa</i> (ID)	49233™
<i>Stenotrophomonas maltophilia</i> (ID)	49130™
MICROSCAN® DRIED GRAM-POSITIVE PANEL QC SET	
<i>Enterococcus faecalis</i> (MIC/ID)	29212™ or 29212-MINI-PACK™
<i>Enterococcus faecalis</i> (MIC)	51299™ or 51299-MINI-PACK™
<i>Enterococcus raffinosus</i> (ID)	49464™
<i>Escherichia coli</i> (MIC)	35218™ or 35218-MINI-PACK™
<i>Micrococcus luteus</i> (ID)	49732™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (MIC/ID/CfxS/lcd)	29213™ or 29213-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (CfxS)	43300™ or 43300-MINI-PACK™
MICROSCAN® DRIED GRAM-POSITIVE PANEL QC SET	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (lcd)	BAA-977™ or BAA-977-MINI-PACK™
<i>Staphylococcus saprophyticus</i> (ID)	49907™
<i>Staphylococcus gallinarum</i> (ID)	49148™
<i>Streptococcus gallolyticus</i> (ID)	49147™
<i>Streptococcus pneumoniae</i> (ID)	49136™
MICROSCAN® DRIED GRAM-NEGATIVE PANEL QC SET	
<i>Enterococcus faecalis</i> (ID)	29212™ or 29212-MINI-PACK™
<i>Escherichia coli</i> (MIC/ID)	25922™ or 25922-MINI-PACK™
<i>Escherichia coli</i> (MIC)	35218™ or 35218-MINI-PACK™
<i>Klebsiella oxytoca</i> (ID)	49131™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> (ESBL)	700603™ or 700603-MINI-PACK™
<i>Proteus vulgaris</i> (ID)	49132™
<i>Pseudomonas aeruginosa</i> (MIC/ID)	27853™ or 27853-MINI-PACK™
<i>Pseudomonas putida</i> (ID)	49128™
<i>Providencia stuartii</i> (ID)	49809™
<i>Ralstonia insidiosa</i> (ID)	49129™
<i>Shewanella haliotis</i> (ID)	49138™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (ID)	29213™ or 29213-MINI-PACK™
MICROSCAN® MICROSTREP PLUS® QC	
<i>Streptococcus pneumoniae</i> (MIC)	49619™ or 49619-MINI-PACK™
MICROSCAN® MICROSTREP PLUS® QC (NON-FDA AND JAPANESE MARKETS ONLY)	
<i>Enterococcus faecalis</i> (ID)	29212™ or 29212-MINI-PACK™
<i>Streptococcus pneumoniae</i> (MIC)	49619™ or 49619-MINI-PACK™

TABLE 2. Beckman Coulter

Sets and Organisms	ATCC® No.
MICROSCAN® RAPID ANAEROBE ID QC SET	
<i>Bacteroides fragilis</i> (ID)	25285™
<i>Clostridium perfringens</i> (ID)	13124™ or 13124-MINI-PACK™
<i>Clostridium sordellii</i> (ID)	9714™ or 9714-MINI-PACK™
<i>Fingoldia magna</i> (ID)	29328™
MICROSCAN® RAPID GRAM-NEGATIVE ID QC SET	
<i>Acinetobacter</i> sp. (ID)	49139™
<i>Aeromonas hydrophila</i> (ID)	49140™
<i>Escherichia coli</i> (ID)	25922™ or 25922-MINI-PACK™
<i>Klebsiella aerogenes</i> (ID)	51697™
<i>Klebsiella oxytoca</i> (ID)	49131™
<i>Shewanella haliotis</i> (ID)	49138™
MICROSCAN® RAPID GRAM-NEGATIVE ID TYPE 4 QC SET	
<i>Klebsiella aerogenes</i> (ID)	51697™
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i> (ID)	49141™
<i>Enterococcus faecalis</i> (ID)	29212™ or 29212-MINI-PACK™
<i>Escherichia coli</i> (ID)	25922™ or 25922-MINI-PACK™
<i>Klebsiella oxytoca</i> (ID)	49131™
<i>Proteus vulgaris</i> (ID)	49132™
<i>Pseudomonas aeruginosa</i> (ID)	27853™ or 27853-MINI-PACK™
<i>Shewanella haliotis</i> (ID)	49138™
MICROSCAN® RAPID GRAM-POSITIVE ID TYPE 2 QC SET	
<i>Acinetobacter</i> sp. (ID)	49139™
<i>Enterococcus casseliflavus</i> (ID)	700327™ or 700327-MINI-PACK™
<i>Enterococcus durans</i> (ID)	6056™
<i>Klebsiella oxytoca</i> (ID)	49131™
<i>Shewanella haliotis</i> (ID)	49138™
<i>Staphylococcus epidermidis</i> (ID)	49134™
<i>Streptococcus pasteurianus</i> (ID)	49133™
MICROSCAN® STREAMLINED ID QC	
<i>Acinetobacter</i> sp. (Rapid Pos ID2)	49139™
<i>Klebsiella aerogenes</i> (Rapid Neg ID4)	51697™
<i>Enterococcus casseliflavus</i> (Rapid Pos ID2)	700327™ or 700327-MINI-PACK™
<i>Enterococcus durans</i> (Rapid Pos ID2)	6056™
<i>Enterococcus faecalis</i> (Rapid Neg ID4)	29212™ or 29212-MINI-PACK™
<i>Escherichia coli</i> (Dried Gram-Negative)	25922™ or 25922-MINI-PACK™
<i>Pseudomonas aeruginosa</i> (Rapid Neg ID4)	27853™ or 27853-MINI-PACK™
<i>Streptococcus gallolyticus</i> (Dried Gram-Positive)	49147™

TABLE 2. Beckman Coulter

Sets and Organisms	ATCC® No.
MICROSCAN® YEAST ID QC SET	
<i>Candida albicans</i> (ID)	66027™
<i>Candida glabrata</i> (ID)	66032™
<i>Candida kefyr</i> (ID)	66028™
<i>Candida lusitanae</i> (ID)	66035™
<i>Candida tropicalis</i> (ID)	66029™
<i>Cryptococcus albidus</i> (ID)	66030™
<i>Papiliotrema laurentii</i> (ID)	66036™
<i>Cryptococcus neoformans</i> (ID)	66031™
<i>Cryptococcus uniguttulatus</i> (ID)	66033™
<i>Rhodotorula mucilaginosa</i> (ID)	66034™

TABLE 3. BiOLOG

Sets and Organisms	ATCC® No.
FF MICROPLATE™ QC SET	
<i>Aspergillus niger</i>	16888™
<i>Galactomyces candidus</i>	34614™
<i>Penicillium chrysogenum</i>	10106™
GEN III MICROPLATE™ QC/ OMNILOG® & MICROSTATION™ VALIDATION SET	
<i>Escherichia coli</i>	11775™ or 11775-MINI-PACK™
<i>Paenibacillus polymyxa</i>	842™
<i>Staphylococcus epidermidis</i>	12228™ or 12228-MINI-PACK™
<i>Stenotrophomonas maltophilia</i>	13637™
AN MICROPLATE™ QC SET	
<i>Bifidobacterium breve</i>	15700™
<i>Clostridium sordellii</i>	9714™ or 9714-MINI-PACK™
<i>Lactobacillus casei</i>	393™
<i>Parvimonas micra</i>	33270™
YT MICROPLATE™ QC SET	
<i>Candida albicans</i>	10231™ or 10231-MINI-PACK™
<i>Candida geochares</i>	36852™
<i>Candida kefyr</i>	2512™
<i>Galactomyces candidus</i>	34614™

TABLE 4. bioMérieux

Sets and Organisms	ATCC® No.
API® 10 S QC SET	
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i>	13047™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Stenotrophomonas maltophilia</i>	51331™
<i>Proteus mirabilis</i>	35659™ or 35659-MINI-PACK™
API® 20A QC SET	
<i>Bacteroides ovatus</i>	8483™
<i>Clostridium perfringens</i>	13124™ or 13124-MINI-PACK™
<i>Clostridium sordellii</i>	9714™ or 9714-MINI-PACK™

TABLE 4. bioMérieux

Sets and Organisms	ATCC® No.
API® 20C AUX QC SET	
<i>Candida glabrata</i>	15126™
<i>Meyerozyma guilliemondii</i>	6260™
<i>Papiliotrema laurentii</i>	18803™
API® 20 E QC SET	
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i>	13047™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	35657™
<i>Proteus mirabilis</i>	35659™ or 35659-MINI-PACK™
<i>Stenotrophomonas maltophilia</i>	51331™
API® 20 NE QC SET	
<i>Aeromonas hydrophila</i>	35654™ or 35654-MINI-PACK™
<i>Alcaligenes faecalis</i> subsp. <i>faecalis</i>	35655™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
<i>Sphingobacterium multivorum</i>	35656™
API® 20 STREP STREAMLINED QC SET	
<i>Streptococcus equi</i> subsp. <i>zooepidemicus</i>	700400™
API® 20 STREP QC SET	
<i>Streptococcus equi</i> subsp. <i>zooepidemicus</i>	700400™
<i>Streptococcus uberis</i>	700407™
API® CAMPY QC SET	
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	49943™
<i>Campylobacter fetus</i> subsp. <i>fetus</i>	25936™
API® CORYNE QC SET	
<i>Cellulosimicrobium cellulans</i>	27402™
<i>Corynebacterium renale</i>	19412™
<i>Microbacterium testaceum</i>	15829™ or 15829-MINI-PACK™
API® RAPID ID 32 A QC SET	
<i>Actinomyces viscosus</i>	15987™
<i>Bacteroides fragilis</i>	23745™
<i>Clostridium sporogenes</i>	19404™ or 19404-MINI-PACK™
<i>Clostridium sordellii</i>	9714™ or 9714-MINI-PACK™
<i>Capnocytophaga sputigena</i>	33612™
API® ID 32 E QC SET	
<i>Vanrja humicola</i>	64676™
<i>Candida glabrata</i>	64677™
API® ID 32 C QC SET	
<i>Escherichia coli</i>	11775™ or 11775-MINI-PACK™
<i>Salmonella enterica</i>	35664™
<i>Raoultella ornithinolytica</i>	31898™
<i>Stenotrophomonas maltophilia</i>	51331™
API® ID 32 GN QC SET	
<i>Delftia acidovorans</i>	43868™
<i>Sphingobacterium multivorum</i>	35656™

TABLE 4. bioMérieux

Sets and Organisms	ATCC® No.
API® ID 32 STAPH QC SET	
<i>Staphylococcus lugdunensis</i>	49576™
<i>Staphylococcus gallinarum</i>	700401™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29213™ or 29213-MINI-PACK™
API® ID 32 STREP QC SET	
<i>Streptococcus equi</i> subsp. <i>equi</i>	33398™
<i>Streptococcus vestibularis</i>	49124™
<i>Streptococcus agalactiae</i>	12401™
API® NH STREAMLINED QC SET	
<i>Neisseria gonorrhoeae</i>	31426™
API® NH QC SET	
<i>Aggregatibacter aphrophilus</i>	49917™
<i>Haemophilus influenzae</i>	10211™
<i>Neisseria gonorrhoeae</i>	31426™
API® STAPH STREAMLINED QC SET	
<i>Staphylococcus capitis</i> subsp. <i>capitis</i>	35661™
API® STAPH QC SET	
<i>Staphylococcus capitis</i> subsp. <i>capitis</i>	35661™
<i>Staphylococcus lentus</i>	700403™
<i>Staphylococcus xylosum</i>	700404™
API® RAPID 20E QC SET	
<i>Escherichia coli</i>	11775™ or 11775-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	35657™
<i>Proteus hauseri</i>	13315™
API® RAPID ID 32 E QC SET	
<i>Klebsiella oxytoca</i>	43863™
<i>Morganella morganii</i> subsp. <i>morganii</i>	25829™
<i>Proteus mirabilis</i>	35659™ or 35659-MINI-PACK™
ETEST® ESBL PM/PML STRIP	
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	700603™ or 700603-MINI-PACK™
ETEST® ESBL CT/CTL STRIP	
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	700603™ or 700603-MINI-PACK™
<i>Escherichia coli</i>	35218™ or 35218-MINI-PACK™
ETEST® ESBL TZ/TZL STRIP	
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	700603™ or 700603-MINI-PACK™
<i>Escherichia coli</i>	35218™ or 35218-MINI-PACK™
ETEST® MBL IP/IPI STRIP	
<i>Stenotrophomonas maltophilia</i>	13636™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
ETEST® MBL MP/MPI STRIP	
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	700603™ or 700603-MINI-PACK™
<i>Klebsiella pneumoniae</i>	BAA-2146™
RAPIDEC STAPH QC SET	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	25923™ or 25923-MINI-PACK™
<i>Staphylococcus saprophyticus</i>	43867™

TABLE 4. bioMérieux

Sets and Organisms	ATCC® No.
VITEK® 2 CBC STREAMLINED QC SET	
<i>Corynebacterium urealyticum</i>	43044™ or 43044-MINI-PACK™
<i>Microbacterium testaceum</i>	15829™ or 15829-MINI-PACK™
VITEK® 2 CBC COMPREHENSIVE QC SET	
<i>Arcanobacterium haemolyticum</i>	BAA-1784™ or BAA-1784-MINI-PACK™
<i>Cellulosimicrobium cellulans</i>	BAA-1816™ or BAA-1816-MINI-PACK™
<i>Cellulosimicrobium cellulans</i>	BAA-1817™ or BAA-1817-MINI-PACK™
<i>Corynebacterium renale</i>	BAA-1785™ or BAA-1785-MINI-PACK™
<i>Corynebacterium urealyticum</i>	43044™ or 43044-MINI-PACK™
<i>Curtobacterium pusillum</i>	19096™ or 19096-MINI-PACK™
<i>Klebsiella oxytoca</i>	700324™ or 700324-MINI-PACK™
<i>Microbacterium liquefaciens</i>	BAA-1819™ or BAA-1819-MINI-PACK™
<i>Microbacterium paraoxydans</i>	BAA-1818™ or BAA-1818-MINI-PACK™
<i>Microbacterium testaceum</i>	15829™ or 15829-MINI-PACK™
<i>Ochrobactrum anthropi</i>	BAA-749™ or BAA-749-MINI-PACK™
VITEK® 2 BCL STREAMLINED QC SET	
<i>Brevibacillus agri</i>	51663™ or 51663-MINI-PACK™
VITEK® 2 BCL COMPREHENSIVE QC SET	
<i>Aneurinibacillus aneurinilyticus</i>	11376™ or 11376-MINI-PACK™
<i>Bacillus badius</i>	14574™ or 14574-MINI-PACK™
<i>Bacillus circulans</i>	61™
<i>Bacillus megaterium</i>	14581™ or 14581-MINI-PACK™
<i>Bacillus pumilus</i>	BAA-1434™ or BAA-1434-MINI-PACK™
<i>Brevibacillus agri</i>	51663™ or 51663-MINI-PACK™
<i>Brevibacillus laterosporus</i>	64™ or 64-MINI-PACK™
<i>Klebsiella aerogenes</i>	13048™ or 13048-MINI-PACK™
<i>Paenibacillus macerans</i>	8509™ or 8509-MINI-PACK™
<i>Paenibacillus polymyxa</i>	7070™ or 7070-MINI-PACK™
<i>Paenibacillus gordonae</i>	29948™ or 29948-MINI-PACK™
<i>Staphylococcus epidermidis</i>	12228™ or 12228-MINI-PACK™
VITEK® 2 GN STREAMLINED QC SET	
<i>Enterobacter hormaechei</i>	700323™ or 700323-MINI-PACK™
<i>Stenotrophomonas maltophilia</i>	17666™
VITEK® 2 GN COMPREHENSIVE QC SET	
<i>Acinetobacter baumannii</i>	BAA-747™ or BAA-747-MINI-PACK™
<i>Enterobacter hormaechei</i>	700323™ or 700323-MINI-PACK™
<i>Elizabethkingia meningoseptica</i>	13253™ or 13253-MINI-PACK™
<i>Klebsiella oxytoca</i>	700324™ or 700324-MINI-PACK™
<i>Ochrobactrum anthropi</i>	BAA-749™ or BAA-749-MINI-PACK™
<i>Proteus vulgaris</i>	6380™ or 6380-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	9721™
<i>Pseudomonas aeruginosa</i>	BAA-1744™ or BAA-1744-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Stenotrophomonas maltophilia</i>	17666™
VITEK® 2 GP STREAMLINED QC SET	
<i>Enterococcus casseliflavus</i>	700327™ or 700327-MINI-PACK™
<i>Staphylococcus saprophyticus</i>	BAA-750™ or BAA-750-MINI-PACK™

TABLE 4. bioMérieux

Sets and Organisms	ATCC® No.
VITEK® 2 GP COMPREHENSIVE QC SET	
<i>Enterococcus casseliflavus</i>	700327™ or 700327-MINI-PACK™
<i>Enterococcus saccharolyticus</i>	43076™ or 43076-MINI-PACK™
<i>Kocuria kristinae</i>	BAA-752™ or BAA-752-MINI-PACK™
<i>Listeria monocytogenes</i>	BAA-751™ or BAA-751-MINI-PACK™
<i>Staphylococcus saprophyticus</i>	BAA-750™ or BAA-750-MINI-PACK™
<i>Staphylococcus sciuri</i> subsp. <i>sciuri</i>	29061™
<i>Streptococcus equi</i> subsp. <i>zooepidemicus</i>	43079™
<i>Streptococcus pneumoniae</i>	49619™ or 49619-MINI-PACK™
<i>Streptococcus thermophilus</i>	19258™ or 19258-MINI-PACK™
VITEK® 2 YST STREAMLINED QC SET	
<i>Candida albicans</i>	14053™ or 14053-MINI-PACK™
VITEK® 2 YST COMPREHENSIVE QC SET	
<i>Candida albicans</i>	14053™ or 14053-MINI-PACK™
<i>Candida glabrata</i>	MYA-2950™ or MYA-2950-MINI-PACK™
<i>Candida lusitanae</i>	34449™ or 34449-MINI-PACK™
<i>Candida utilis</i>	9950™ or 9950-MINI-PACK™
<i>Hanseniaspora valbyensis</i>	58370™
<i>Oligella ureolytica</i>	43534™ or 43534-MINI-PACK™
<i>Prototheca wickerhamii</i>	16529™
<i>Sporidiobolus salmonicolor</i>	MYA-4550™
<i>Staphylococcus epidermidis</i>	12228™ or 12228-MINI-PACK™
<i>Cutaneotrichosporon dermatis</i>	204094™
<i>Zygosaccharomyces parabailii</i>	MYA-4549™ or MYA-4549-MINI-PACK™
VITEK® 2 NH STREAMLINED QC SET	
<i>Eikenella corrodens</i>	BAA-1152™ or BAA-1152-MINI-PACK™
VITEK® 2 NH COMPREHENSIVE QC SET	
<i>Aggregatibacter aphrophilus</i>	33389™ or 33389-MINI-PACK™
<i>Eikenella corrodens</i>	BAA-1152™ or BAA-1152-MINI-PACK™
<i>Klebsiella aerogenes</i>	13048™ or 13048-MINI-PACK™
<i>Haemophilus influenzae</i>	9007™ or 9007-MINI-PACK™
<i>Neisseria gonorrhoeae</i>	19424™ or 19424-MINI-PACK™
<i>Neisseria lactamica</i>	23970™ or 23970-MINI-PACK™
<i>Oligella urethralis</i>	17960™ or 17960-MINI-PACK™
<i>Paenibacillus polymyxa</i>	7070™ or 7070-MINI-PACK™
<i>Staphylococcus epidermidis</i>	12228™ or 12228-MINI-PACK™
VITEK® 2 ANC STREAMLINED QC SET	
<i>Bacteroides ovatus</i>	BAA-1296™ or BAA-1296-MINI-PACK™
<i>Clostridium septicum</i>	12464™ or 12464-MINI-PACK™
VITEK® 2 ANC COMPREHENSIVE QC SET	
<i>Bacteroides ovatus</i>	BAA-1296™ or BAA-1296-MINI-PACK™
<i>Bacteroides vulgatus</i>	8482™ or 8482-MINI-PACK™
<i>Clostridium perfringens</i>	13124™ or 13124-MINI-PACK™
<i>Clostridium septicum</i>	12464™ or 12464-MINI-PACK™
<i>Clostridium sordellii</i>	9714™ or 9714-MINI-PACK™
<i>Corynebacterium striatum</i>	BAA-1293™ or BAA-1293-MINI-PACK™
<i>Parabacteroides distasonis</i>	BAA-1295™ or BAA-1295-MINI-PACK™

TABLE 4. bioMérieux

Sets and Organisms	ATCC® No.
VITEK® 2 AST-GN QC SET	
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Escherichia coli</i>	35218™ or 35218-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	700603™ or 700603-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
VITEK® 2 AST-GP QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Enterococcus faecalis</i>	51299™ or 51299-MINI-PACK™
<i>Escherichia coli</i>	35218™ or 35218-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29213™ or 29213-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	BAA-976™ or BAA-976-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	BAA-977™ or BAA-977-MINI-PACK™
<i>Staphylococcus aureus</i>	BAA-1026™ or BAA-1026-MINI-PACK™
<i>Streptococcus pneumoniae</i>	49619™ or 49619-MINI-PACK™
VITEK® 2 AST-YS QC SET	
<i>Issatchenkia orientalis</i>	6258™ or 6258-MINI-PACK™*
<i>Candida parapsilosis</i>	22019™ or 22019-MINI-PACK™

TABLE 5. CHARM® Sciences

Sets and Organisms	ATCC® No.
E*COLITE™ TEST/PATHOGEL™ TEST	
<i>Enterobacter cloacae</i> subsp. <i>cloacae</i>	35030™
<i>Escherichia coli</i>	11775™ or 11775-MINI-PACK™
<i>Klebsiella variicola</i>	31488™ or 31488-MINI-PACK™
<i>Stenotrophomonas maltophilia</i>	13637™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium (Pathogel only)	14028™ or 14028-MINI-PACK™
FAST PHAGE™ TEST	
<i>Escherichia coli</i>	700609™
<i>Escherichia coli</i>	700891™
<i>Escherichia coli</i> bacteriophage Phi X174	13706-B1™
<i>Escherichia coli</i> bacteriophage MS2	15597-B1™

TABLE 6. CPI International

Sets and Organisms	ATCC® No.
COLITAG™ TEST KIT QC SET	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium	700720™
<i>Pseudomonas aeruginosa</i>	10145™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	13883™
<i>Citrobacter freundii</i>	8090™
<i>Escherichia coli</i>	11775™ or 11775-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™

TABLE 7. IDEXX

Sets and Organisms	ATCC® No.
COLILERT® QC SET	
<i>Escherichia coli</i>	25922™, 11775™, 25922-MINI-PACK™, or 11775-MINI-PACK™
<i>Klebsiella variicola</i>	31488™ or 31488-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	10145™, 27853™, or 27853-MINI-PACK™
COLILERT® 18 QC SET	
<i>Escherichia coli</i>	11775™ or 11775-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Klebsiella variicola</i>	31488™ or 31488-MINI-PACK™
<i>Klebsiella aerogenes</i>	13048™ or 13048-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	10145™, 27853™, or 27853-MINI-PACK™
COLILERT® 250 QC SET	
<i>Escherichia coli</i>	25922™, 11775™, 25922-MINI-PACK™, or 11775-MINI-PACK™
<i>Klebsiella variicola</i>	31488™ or 31488-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	10145™, 27853™, or 27853-MINI-PACK™
COLISURE® QC SET	
<i>Escherichia coli</i>	25922™, 11775™, 25922-MINI-PACK™, or 11775-MINI-PACK™
<i>Klebsiella variicola</i>	31488™ or 31488-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	10145™, 27853™, or 27853-MINI-PACK™
ENTEROLERT® QC SET	
<i>Enterococcus faecium</i>	35667™
<i>Serratia marcescens</i>	43862™
<i>Aerococcus viridans</i>	10400™
ENTEROLERT® E QC SET	
<i>Aerococcus viridans</i>	10400™
<i>Enterococcus faecium</i>	35667™
<i>Serratia marcescens</i>	43862™
ENTEROLERT® DW QC SET	
<i>Enterococcus faecium</i>	35667™
<i>Serratia marcescens</i>	43862™
ENTEROLERT® 250 QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Aerococcus viridans</i>	11563™
LEGIOLERT® QC SET	
<i>Legionella pneumophila</i> subsp. <i>pneumophila</i>	33152™
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
PSEUDALERT® QC SET	
<i>Pseudomonas aeruginosa</i>	27853™, 27853-MINI-PACK™, or 10145™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Pseudomonas fluorescens</i>	13525™
PSEUDALERT® 250 QC SET	
<i>Pseudomonas aeruginosa</i>	27853™, 27853-MINI-PACK™, or 10145™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Pseudomonas fluorescens</i>	13525™


TABLE 8. Thermo Fisher Scientific


Sets and Organisms	ATCC® No.
RapID™ ANA II PANEL QC Set	
<i>Bacteroides uniformis</i>	8492™
<i>Clostridium sordellii</i>	9714™ or 9714-MINI-PACK™
<i>Parabacteroides distasonis</i>	8503™
RapID™ CB PLUS QC SET	
<i>Trueperella pyogenes</i>	19411™
<i>Corynebacterium pseudodiphtheriticum</i>	10701™
<i>Paenibacillus polymyxa</i>	842™
RapID™ NH QC SET	
<i>Haemophilus influenzae</i> Type a	9006™
<i>Aggregatibacter aphrophilus</i>	49146™
<i>Moraxella catarrhalis</i>	8176™
<i>Oligella urethralis</i>	17960™ or 17960-MINI-PACK™
RapID™ NF PLUS QC SET	
<i>Acinetobacter baumannii</i>	19606™ or 19606-MINI-PACK™
<i>Aeromonas hydrophila</i>	35654™ or 35654-MINI-PACK™
<i>Elizabethkingia meningoseptica</i>	13253™ or 13253-MINI-PACK™
<i>Oligella ureolytica</i>	43534™ or 43534-MINI-PACK™
RapID™ ONE QC SET	
<i>Klebsiella aerogenes</i>	13048™ or 13048-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Proteus vulgaris</i>	6380™ or 6380-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
RapID™ SS/U QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	13883™
<i>Proteus mirabilis</i>	25933™
<i>Serratia marcescens</i>	8100™
RapID™ STAPH PLUS QC SET	
<i>Staphylococcus haemolyticus</i>	29970™
<i>Staphylococcus saprophyticus</i>	35552™
<i>Klebsiella aerogenes</i>	13048™ or 13048-MINI-PACK™
<i>Oligella ureolytica</i>	43534™ or 43534-MINI-PACK™
RapID™ STR QC SET	
<i>Enterococcus durans</i>	11576™
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Streptococcus gallolyticus</i>	9809™
<i>Streptococcus pyogenes</i>	19615™ or 19615-MINI-PACK™
RapID™ YEAST PLUS QC SET	
<i>Candida albicans</i>	14053™ or 14053-MINI-PACK™
<i>Candida glabrata</i>	2001™
<i>Candida kefyr</i>	2512™
<i>Papiliotrema laurentii</i>	66036™
<i>Yarrowia lipolytica</i>	9773™

TABLE 8. Thermo Fisher Scientific


Sets and Organisms	ATCC® No.
SENSITITRE™ ANAEROBE MIC SUSCEPTIBILITY QC SET	
<i>Bacteroides fragilis</i>	25285™
<i>Bacteroides thetaiotaomicron</i>	29741™
<i>Eggerthella lenta</i>	43055™
SENSITITRE™ AQUACULTURE MIC QC SET	
<i>Aeromonas salmonicida</i> subsp. <i>salmonicida</i>	33658™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
SENSITITRE™ MIC BREAKPOINT QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Escherichia coli</i>	35218™ or 35218-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29213™ or 29213-MINI-PACK™
SENSITITRE™ AUTOMATED GRAM-NEGATIVE ID QC SET	
<i>Edwardsiella tarda</i>	15947™
<i>Klebsiella oxytoca</i>	8724™
<i>Morganella morganii</i> subsp. <i>morganii</i>	25830™
<i>Proteus vulgaris</i>	6896™
<i>Pseudomonas aeruginosa</i>	10145™
<i>Shigella sonnei</i>	25931™ or 25931-MINI-PACK™
SENSITITRE™ AUTOMATED GRAM-POSITIVE ID QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Kocuria rosea</i>	186™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
<i>Staphylococcus epidermidis</i>	700296™
SENSITITRE™ H. INFLUENZAE/S. PNEUMONIAE MIC SUSCEPTIBILITY QC SET	
<i>Haemophilus influenzae</i>	49247™
<i>Haemophilus influenzae</i>	49766™
<i>Streptococcus pneumoniae</i>	49619™ or 49619-MINI-PACK™
SENSITITRE™ MIC SUSCEPTIBILITY QC SET	
<i>Enterococcus faecalis</i>	29212™ or 29212-MINI-PACK™
<i>Enterococcus faecalis</i>	51299™ or 51299-MINI-PACK™
<i>Escherichia coli</i>	35218™ or 35218-MINI-PACK™
<i>Escherichia coli</i>	25922™ or 25922-MINI-PACK™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	700603™ or 700603-MINI-PACK™
<i>Pseudomonas aeruginosa</i>	27853™ or 27853-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	29213™ or 29213-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	BAA-976™ or BAA-976-MINI-PACK™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	BAA-977™ or BAA-977-MINI-PACK™
SENSITITRE™ YEAST ONE QC SET	
<i>Candida parapsilosis</i>	22019™ or 22019-MINI-PACK™
<i>Issatchenkia orientalis</i>	6258™ or 6258-MINI-PACK™

CREDIBLE STANDARDS INCREDIBLE RESULTS

 10801 University Boulevard
Manassas, Virginia 20110-2209

 703.263.8277

 703.365.2701

 sales@atcc.org

 www.atcc.org

QCS-10092019-16

©2019 American Type Culture Collection. The ATCC trademark and trade name, and any other trademarks listed in this publication are trademarks owned by the American Type Culture Collection unless indicated otherwise. BD Crystal™, Enterotube™, Phoenix™, and GeneOhm™ are trademarks of Becton, Dickinson and Co. FF Microplate™ and Microstation™ are trademarks of BiOLOG, Inc. Omnilog® is a registered trademark of BiOLOG, Inc. API®, Etest® and VITEK® are registered trademarks of BioMérieux, Inc. E*Colite™, Fast Phage™, and Pathogel™ are trademarks of Charm Sciences, Inc. Colitag™ is a trademark of CPI International. Colilert®, Colisure® and Enterolert® are registered trademarks of IDEXX Laboratories, Inc. RapID™ and Sensititre™ are trademarks of Thermo Fisher Scientific. Beckman Coulter, MicroScan® and MICroSTREP plus® are trademarks of Beckman Coulter, Inc.

These products are for laboratory use only. Not for human or diagnostic use. ATCC products may not be resold, modified for resale, used to provide commercial services or to manufacture commercial products without prior ATCC written approval.