

ATCC Medium: 2688 ATCC medium: Acidovorax Complex Medium (ACM)

KH ₂ PO ₄ ,	0.30 g
MgSO ₄	0.30 g
MOPS.....	2.00 g
Yeast extract.....	0.50 g
Pyruvic acid,.....	2.00 g
L-glutamine.....	2.00 g
Distilled water.....	1.0 L

Adjust the pH to 7.2 to 7.5. Autoclave at 121C for 15 minutes

For solid media, use 1.5% Phytigel (SigmaAldrich agar or agarose will not support growth)