

ATCC medium: 2720 Carboxymethyl Cellulose Medium

(NH ₄) ₂ SO ₄	1.0g	
MgSO ₄ x 7H ₂ O	1.0g	
CaCl ₂ x 2H ₂ O	1.0g	
FeCl ₃	0.2g	
K ₂ HPO ₄ (Filter)	1.0g	
Casitone	2.0g	
Carboxymethyl Cellulose	15.0g	
Agar	15.0g	
Distilled water		950 ml

Autoclave, then aseptically add K₂HPO₄ 50ml solution.
Final pH is 7.0