

ATCC medium: 187 *Pseudomonas* medium I

Tryptone (BD 211705).....10.0 g
Yeast extract.....10.0 g
K₂HPO₄5.0 g
Glucose.....5.0 g
Salts Solution (see below)...5.0 ml
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.

Salts Solution:

MgSO₄ . 7H₂O4.0 g
NaCl.....0.2 g
FeSO₄0.2 g
MnSO₄ . 4H₂O0.2 g
Distilled water.....100.0 ml