

**ATCC medium: 1394 *Halobacterium pharaonis* medium**

Casamino Acids (BD 223050).....15.0 g  
Sodium citrate . 2H<sub>2</sub>O.....3.0 g  
Glutamic acid.....2.5 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....2.5 g  
KCl.....2.0 g  
NaCl.....250.0 g  
Agar.....20.0 g  
Distilled water.....1.0 L

Adjust pH prior to sterilization to 6.0. Autoclave at 121C for 15 minutes. Cool to 50C. Re-adjust pH aseptically to 8.5 after autoclaving.