

ATCC Medium 129: Nutrient Agar/Stabs/Broth with 0.5% NaCl

Nutrient Agar (BD 213000)..... 23.0 g
NaCl (JT Baker 3624.01) 5.0 g
DI Water..... 1000.0 mL

Autoclave at 121°C.

For Stabs

Nutrient Broth (BD 234000)..... 8.0 g
NaCl (JT Baker 3624.01)..... 5.0 g
Agar (BD 214010)..... 5.0 g
DI Water..... 1000.0 mL

Autoclave at 121°C.

For Broth

Nutrient Broth (BD 234000)..... 8.0 g
NaCl (JT Baker 3624.01)..... 5.0 g
DI Water..... 1000.0 mL

Autoclave at 121°C.