

ATCC Medium: 2343 SME-SAN2 medium

NaCl	27.7 g
MgSO ₄ . 7H ₂ O.....	7.0 g
MgCl ₂ . 6H ₂ O.....	5.5 g
KCl.....	0.65 g
NaBr.....	100.0 mg
H ₃ BO ₃	30.0 mg
SrCl . 6H ₂ O.....	15.0 mg
KI.....	0.05 mg
CaCl ₂ . 2H ₂ O.....	1.5 g
Yeast extract.....	1.0 g
BisTris/HCl pH 7.0 (1 M solution)	25.0 ml
Distilled H ₂ O to	1.0 L

Autoclave at 121C for 15 minutes.