

ATCC Medium 252: Brain Heart Infusion with 1.2% Agar, 0.5% Yeast Extract, 10% Sucrose, 10% Inactivated horse Serum and 100units/ml Penicillin

Base Medium

Brain Heart Infusion	37.0 g
Agar (if required).....	12.0 g
Yeast Extract	5.0 g
Sucrose	100.0 g
DI Water.....	850.0 mL

Autoclave at 121°C. Cool to ~ 50-55°C

Aseptically add the following two components:

Horse Serum (heat inactivated).....	100.0 mL
-------------------------------------	----------

Penicillin solution (0.063 mg/mL)

Penicillin (Sigma, P7794)	0.063 g
DI Water.....	50.0 mL

Filter sterilize.