

**ATCC medium: 1235 LB Medium (ATCC medium 1065) with 5 mcg/ml tetracycline and 10 mcg/ml ampicillin**

ATCC Medium 1065 (see below) with 5 mcg/ml tetracycline and 10 mcg/ml ampicillin

*ATCC Medium 1065:*

Tryptone (BD 211705).....	10.0 g
Yeast Extract (BD 212750).....	5.0 g
NaCl.....	10.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.