

ATCC medium: 2136 CY medium

Casitone.....5.0 g
Yeast extract.....1.0 g
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.

NOTE: Growth of *Zoogloea resiniphila* is poor on this medium, however it produces a dispersed turbid suspension that can be freeze-dried.

A better growth medium would be tryptic soy broth (30.0 g/L) with addition of filter-sterilized sodium pyruvate (1.0 g/L). Incubate at 30C in a roller tube. Growth is in the form of a large flake or clump due to copious production of exopolymer which interferes with lyophilization.