

**ATCC medium: 840 Wilbrinck's agar for *Xanthomonas albilineans***

Sucrose.....10.0 g  
Peptone.....5.0 g  
K<sub>2</sub>HPO<sub>4</sub> .....0.5 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....0.25 g  
Na<sub>2</sub>SO<sub>3</sub> (anhydrous) .....0.05 g  
Agar.....20.0 g  
Distilled water.....1.0 L

Adjust final pH to 7.2.