

**ATCC Medium: 766 *Desulfovibrio* Medium with 1% Lactate**

Peptone.....	5.0 g
Beef Extract.....	3.0 g
Yeast Extract.....	0.2 g
MgSO <sub>4</sub> .....	1.5 g
Na <sub>2</sub> SO <sub>4</sub> .....	1.5 g
Fe(NH <sub>4</sub> ) <sub>2</sub> (SO <sub>4</sub> ) .....	0.1 g
Glucose.....	5.0 g
Lactate.....	10.0 g
Agar.....	15.0 g
Tap Water.....	1000 ml

Adjust pH to 7.0 +/- 0.2. Autoclave medium at 121°C, cool under N<sub>2</sub> and dispense anaerobically.