




Product Sheet

Pneumococcal

Polysaccharide Powder

Type 3 (US Type 3) (ATCC® 32-X™)

Please read this FIRST



Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: *Pneumococcal Polysaccharide Powder Type 3 (US Type 3) (ATCC® 32-X™)*

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor

Description

Antigenic Properties: Serotype 3 (Danish designation 3)

Product Description: Purified pneumococcal polysaccharide type 3 is manufactured by Pfizer, Inc. for small quantity distribution by ATCC. The pneumococcal powder type 3 is a component of the current 13-valent vaccine against *Streptococcus pneumoniae*.

The pneumococcal polysaccharide type 3 can be used to assay for type-specific antibodies, research immune responses and study immunodeficiencies.

Propagation

Method for Resuspension

1. Prepare a 2 mg/mL pneumococcal polysaccharide stock solution by adding 1 mL of pyrogen-free distilled water to 2 mg of pneumococcal polysaccharide.
2. Incubate the suspension overnight at 2-8°C with gentle stirring on a stir plate to dissolve the polysaccharide.
3. Aliquot the polysaccharide solution and store the aliquots containing the stock solution at -60°C to -80°C.

Notes

Stock solution can be diluted to a working concentration in the buffer of the researcher's choice.

Storage of Powder: Store at -20°C or colder.

Storage of Reconstituted Material: Store at -60°C to -80°C.

Quantities available: (See Web site for current prices)

2 mg

10 mg

200 mg

Pneumococcal polysaccharide powder type 3 is manufactured by Pfizer, Inc. This is **NOT** a sterile product. Antigenicity is verified by Pfizer, Inc.

References

References and other information relating to this product are available online at www.atcc.org.

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

ATCC Warranty

The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for



Product Sheet

***Pneumococcal
Polysaccharide Powder
Type 3 (US Type 3) (ATCC®
32-X™)***

Please read this **FIRST**



Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: *Pneumococcal Polysaccharide Powder Type 3 (US Type 3) (ATCC® 32-X™)*

damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at www.atcc.org

Additional information on this culture is available on the ATCC web site at www.atcc.org.

© ATCC 2013. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [07/09]

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor