



Product Sheet

HPV-18 purified plasmid DNA (ATCC® 45152D™)

human papilloma virus type 18

Please read this **FIRST**

Storage Temp.
Store at 4°C

Biosafety Level
2

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: *human papilloma virus type 18* (ATCC® 45152D™)

Shipping Information

200 ng purified plasmid DNA, dried. Rehydrate with TE or water.

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor

Description

This is a clone of human papilloma virus type 18 that has been linearized at the unique EcoRI site (nt 2440) and cloned into the EcoRI site of pBR322.

Designation: HPV-18 purified plasmid DNA

Notes

Restriction digests of the clone give the following sizes (kb): BamHI--6.6, 4.8, 1.1; EcoRI--8.0, 4.5; PstI--6.6, 4.4, 0.55, 0.5(doublet).

- ATCC staff

The depositor has requested that this deposit be distributed exclusively for scientific and not commercial purposes.

Insert Information

Insert size (kb): 7.857

Genome: human papilloma virus type 18

Gene name: genome, full-length

Contains complete coding sequence?: Y

Type of DNA: genomic

Insert ends: EcoRI

DNA Seq. Acc.: X05015

Vector Information

Size of Construct (kb): 12.2

Excise insert: EcoRI

Vector: pBR322

Size of intact vector (kb): 4.363

Markers: ampR

Vector ends: EcoRI

References

References and other information relating to this product are available online at www.atcc.org.

Biosafety Level: 2

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

ATCC Warranty

The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While



Product Sheet

HPV-18 purified plasmid DNA (ATCC® 45152D™)

human papilloma virus type 18



reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at www.atcc.org

Additional information on this culture is available on the ATCC web site at www.atcc.org.

© ATCC 2014. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [08/13]

Please read this FIRST

	Storage Temp. Store at 4°C
.....	
	Biosafety Level 2

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: *human papilloma virus type 18* (ATCC® 45152D™)

Shipping Information

200 ng purified plasmid DNA, dried. Rehydrate with TE or water.

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor