



Nucleic Acid Product Sheet

Micrococcus luteus (Schroeter) Cohn (ATCC® 53598D™)

Please read this **FIRST**



Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: *Micrococcus luteus* (Schroeter) Cohn (ATCC® 53598D™)

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor

Description

Source: *Micrococcus luteus* (Schroeter) Cohn

Designation: Genomic DNA from SC 14641

Description:

Genomic DNA isolated from *Micrococcus luteus*. This bacterial strain is also available as ATCC Catalog No.: 53598.

Note:

Genomic DNA isolated from bacteria is appropriate for PCR* and other molecular biology applications.

*The polymerase chain reaction (PCR) process is covered by patents owned by Hoffmann-LaRoche Inc. Use of the PCR process requires a license.

Batch-Specific Information

Refer to the Certificate of Analysis for batch-specific test results.

Preparation Procedure

1. Rehydrate contents of vial with 1X TE or molecular grade H₂O.
2. Place vial at 37°C for 1 hour or at 4°C overnight.

Quality Control Information

QUALITY CONTROL INFORMATION:

1. Bacterial genomic DNA is provided in a dried form. Store at 4°C upon receipt (-20°C if stored for more than 6 months). **Note:** Do not store in freezers with a defrost cycle. Subjecting genomic DNA to repeated freeze/thaw cycles may result in shearing of the DNA.
2. Integrity of DNA was determined by electrophoresis on a 1% agarose gel stained with ethidium bromide, and was found to be intact and of high molecular weight.
3. No RNA was detected by electrophoresis.

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

ATCC Warranty

The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at www.atcc.org

Additional information on this culture is available on the ATCC web site at www.atcc.org.

© ATCC 2012. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [12/05]