



Product Sheet

pRE112 (ATCC® 87692™)

Please read this FIRST



Storage Temp.
**Store frozen vial
at -80⁰ C.**



Biosafety Level
1

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pRE112 (ATCC® 87692™)

Shipping Information

Frozen glycerol stock of *E. coli* containing the plasmid

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor

Description

One of 4 allelic exchange suicide vectors (ATCC# 87691-87694) that provide both selection for chromosomal integration (ampR, cmLR, kanR or tetR) and counterselection for loss of vector DNA and the wild type allele. The construct is a suicide plasmid in any host not expressing *pir*. Cloned inserts may be integrated into the host chromosome (a single recombination event) following electroporation and appropriate antibiotic selection. - Gene 207: 149-157, 1998. The conditional R6K origin of replication requires that the pi protein be expressed in trans for plasmid maintenance. - Cell 15: 1199-1208, 1978. Negative selection for sucrose sensitivity (sacB) selects for a second recombination event resulting in loss of vector DNA. The sacB1 allele is a variant of sacB with certain restriction sites removed by site directed mutagenesis. Expression of sacB confers sensitivity to sucrose. Mol. Microbiol. 18: 877-889, 1995.
Designation: pRE112

Media:

ATCC® Medium 1675: LB medium (ATCC medium 1065) with 10 mcg/ml Chloramphenicol

Conditions

Temperature: 37.0°C

Vector Information

Size (kb): 6.0000000000000000
Vector: pRE112 (plasmid)
Construction: pRE107, cmIR (pACYC184)
Marker(s): cmIR, sacB
Construct size (kb): 6.0
Features: marker(s): cmIR
marker(s): sacB1
replicon: conditional R6K ori
replicon: oriT
MCS: HindIII...EcoRI
restriction site: EcoRI HindIII

References

References and other information relating to this product are available online at www.atcc.org.

Notes

Restriction digests of the clone gave the following sizes (in kb): KpnI 6.55; XbaI 6.55; HindIII 3.2, 3.0.
ATCC Staff

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

ATCC Warranty

The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans. While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has



Product Sheet

pRE112 (ATCC® 87692™)

been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at www.atcc.org

Additional information on this culture is available on the ATCC web site at www.atcc.org.
© ATCC 2013. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [02/22]

Please read this FIRST

Storage Temp.
**Store frozen vial
at -80⁰ C.**

Biosafety Level
1

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pRE112 (ATCC® 87692™)

Shipping Information

Frozen glycerol stock of *E. coli* containing the plasmid

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor