

Antibiotic Susceptibility of *Pseudomonas paraeruginosa* ATCC® BAA-3270™

Antimicrobial	MIC ^a	Interpretation ^b	Antimicrobial	MIC ^a	Interpretation ^b
Amikacin	≥64	R	Gentamicin	≥16	R
Cefazolin	≥64	R	Imipenem	≥16	R
Cefepime	≥64	R	Levofloxacin	≥8	R
Ceftazidime	≥64	R	Meropenem	≥16	R
Ceftazidime / Avibactam	≥16	R	Piperacillin / Tazobactam	≥128	R
Ceftolozane / Tazobactam	≥32	R	Tobramycin	≥16	R
Ciprofloxacin	≥4	R	Ticarcillin / Clavulanic Acid	≥128	R

^a Antibiotic susceptibility was obtained using Vitek2 AST-GN69 and AST-XN09 cards, run as *Pseudomonas aeruginosa*.

^b Parameter set: MIC Interpretation Guideline: CLSI M100-S27 (2017)
Therapeutic Interpretation Guideline: Natural Resistance