

**ATCC medium: 1029 LIT medium**

Liver Infusion Broth (BD 226920).....9.0 g  
Tryptose (BD 211713).....5.0 g  
NaCl.....1.0 g  
Na<sub>2</sub>HPO<sub>4</sub> .....8.0 g  
KCl.....0.4 g  
Glucose.....1.0 g  
Fetal bovine serum (heat-inactivated)..100.0 ml  
Hemin.....10.0 mg  
Distilled water to.....1.0 L

Adjust pH to 7.2 and filter-sterilize. Dispense aseptically in 5.0 ml aliquots into 16 X 125 screw-capped test tubes.