

ATCC Medium 1176: Halobacterium Medium

NaCl	156.0 g
MgCl ₂ x 6H ₂ O	13.0 g
MgSO ₄ x 7H ₂ O	20.0 g
CaCl ₂ x 2H ₂ O	1.0 g
KCl	4.0 g
NaHCO ₃	0.2 g
NaBr	0.5 g
Yeast Extract	5.0 g
Glucose	1.0 g
Agar (if required)	15.0 g
DI Water	1.0 L

Adjust pH to 7.0.

Autoclave at 121°C.