

ATCC medium: 1274 Benzoate medium II

Sodium benzoate.....1.0 g
Yeast extract.....0.5 g
(NH₄)₂HPO₄.....3.0 g
KH₂PO₄1.2 g
NaCl.....1.67 g
FeSO₄ . 7H₂O0.1 g
MgSO₄ . 7H₂O0.2 g
Agar, Noble (BD 214230).....30.0 g
Distilled water.....1.5 L

To retard precipitation of salts, autoclave agar as a separate solution.

Autoclave at 121C for 15 minutes.