

ATCC medium: 1463 Modified Hayflick medium

Heart Infusion Broth (BD 238400).....28.5 g
Agar, Noble (BD 214230), if necessary....9.6 g
Distilled water.....900.0 ml

Autoclave at 121C for 15 minutes. Cool to 50C and add aseptically:

Horse serum (normal).....200.0 ml
Yeast Extract Solutuion (GIBCO 18180)..100.0 ml
0.2% Calf Thymus DNA (Sigma D1501).....12.0 ml

Adjust medium to pH 7.8.