

**ATCC medium: 1487 Marine salts medium for *Sporosarcina halophila***

Proteose Peptone No. 3 (BD 211693).....	5.0 g
Yeast extract.....	10.0 g
Glucose.....	1.0 g
NaCl.....	71.5 g
MgSO <sub>4</sub> . 7H <sub>2</sub> O .....	17.6
MgCl <sub>2</sub> . 6H <sub>2</sub> O .....	14.0 g
CaCl <sub>2</sub> . 2H <sub>2</sub> O .....	3.6 g
KCl.....	1.9 g
NaHCO <sub>3</sub> .....	0.5 g
Agar (if needed).....	20.0 g
Distilled water to.....	1.0 L

Adjust medium for final pH 7.0. Autoclave at 121C for 15 minutes.