

ATCC medium: 1550 Glycocholate mineral medium

K ₂ HPO ₄	3.5 g
KH ₂ PO ₄	1.5 g
(NH ₄) ₂ SO ₄	2.0 g
MgSO ₄ . 7H ₂ O	0.1 g
CaCl ₂ . 2H ₂ O	0.01 g
FeSO ₄ . 7H ₂ O	0.5 mg
Yeast extract.....	0.1 g
Sodium glycocholate.....	2.0 g
Agar.....	15.0 g
Deionized water.....	1.0 L

Adjust pH to 7.0. Autoclave at 121C for 15 minutes.