

ATCC medium: 1663 Dilute potato medium

Potato Decoction (see below).....100.0 ml
Peptone.....0.05 g
Glucose.....1.0 g
Ca(NO₃)₂ · 4H₂O.....0.05 g
Na₂HPO₄0.12 g
Distilled water.....990.0 ml

Adjust medium for final pH 6.8. Autoclave at 121C for 15 minutes.

Potato Decoction:

Diced potato.....20.0 g
Distilled water.....1.0 L

Boil for 30 minutes; filter out potato and return volume of liquid to 1.0 L.