

ATCC medium: 1689 MH medium

Yeast Extract (BD 212750).....	10.0 g
Proteose Peptone No. 3 (BD 211693).....	5.0 g
Glucose.....	1.0 g
NaCl.....	60.7 g
MgCl ₂ . 6H ₂ O	15.0 g
MgSO ₄ . 7H ₂ O	7.4 g
CaCl ₂	0.27 g
KCL.....	1.5 g
NaHCO ₃	45.0 mg
NaBr.....	19.0 mg
Agar.....	20.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.