

ATCC medium: 1708 HM medium

NaCl.....	81.0 g
MgCl ₂	7.0 g
MgSO ₄	9.6 g
CaCl ₂	0.36 g
KCl.....	2.0 g
NaHCO ₃	60.0 mg
NaBr.....	26.0 mg
Proteose Peptone No. 3 (BD 211693).....	5.0 g
Yeast extract.....	10.0 g
Glucose.....	1.0 g
Distilled water to.....	1.0 L

Adjust for final pH 7.1, +/- 0.1. Autoclave at 121C for 15 minutes.