

ATCC Medium: 1724 Revised Ducreyi Medium

Solution A

Heart Infusion Broth (BD 238400).....	25.0 g
Agar (if necessary).....	15.0 g
DI Water.....	400.0 ml

Autoclave at 121°C.

Solution B

Hemoglobin, Dried Bovine (BD 212392).....	10.0 g
DI Water.....	500.0 ml

Autoclave at 121°C.

Solution C

Fetal Bovine Serum.....	100.0 ml
IsoVitalex Enrichment (BD 211876).....	10.0 ml

Filter-sterilize and aseptically combine Solutions A, B, and C at 50°C or less.