

ATCC medium: 1787 KC medium

Broth:

Tryptone.....10.0 g
KCl.....5.0 g
1 M CaCl₂ (add aseptically after autoclaving medium)0.5 ml
Distilled water.....1.0 L

Top Agar:

Tryptone.....10.0 g
NaCl.....5.0 g
Agar.....8.0 g
Distilled water.....1.0 L

Bottom Agar:

Tryptone.....10.0 g
NaCl.....2.5 g
KCl.....2.5 g
1 M CaCl₂ (add aseptically after autoclaving medium)1.0 ml
Agar.....10.0 g
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.