

**ATCC medium: 1826 PD 2 medium**

|                                                            |          |
|------------------------------------------------------------|----------|
| Tryptone (BD 211705).....                                  | 4.0 g    |
| Soytone (BD 243620).....                                   | 2.0 g    |
| Sodium citrate.....                                        | 1.0 g    |
| Sodium succinate.....                                      | 1.0 g    |
| 0.1% Hemin chloride (dissolved in 0.05 N NaOH).....        | 10.0 ml  |
| MgSO <sub>4</sub> . 7H <sub>2</sub> O .....                | 1.0 g    |
| K <sub>2</sub> HPO <sub>4</sub> .....                      | 1.5 g    |
| KH <sub>2</sub> PO <sub>4</sub> .....                      | 1.0 g    |
| Agar.....                                                  | 15.0 g   |
| 20% Aqueous bovine serum albumin fraction V (see below)... | 10.0 ml  |
| Distilled water.....                                       | 980.0 ml |

Dissolve all components except bovine serum albumin (BSA-V) in distilled water and autoclave at 121C for 15 minutes. Cool to approximately 50C and aseptically add filter-sterilized BSA-V solution. Final pH 7.0.