

ATCC medium: 1863 Halophile medium (NCIMB 219)

Casamino Acids (BD 223050).....7.50 g
Yeast extract.....10.0 g
Trisodium citrate.....3.0 g
KCl.....2.0 g
MgSO₄ . 7H₂O20.0 g
FeCl₂ . 4H₂O36.0 mg
MnCl₂ . 4H₂O0.36 mg
NaCl.....250.0 g
Agar (if necessary).....20.0 g
Distilled water.....1.0 L

Adjust medium for final pH 7.4. Autoclave at 121C for 15 minutes.