

ATCC medium: 186 *Pseudomonas* medium

Glucose.....	10.0 g
Yeast extract.....	10.0 g
Peptone.....	1.0 g
K ₂ HPO ₄	1.0 g
KH ₂ PO ₄	0.5 g
NH ₄ Cl.....	0.5 g
Na ₂ SO ₃	0.1 g
CaCl ₂	0.1 g
NaHCO ₃	0.1 g
FeCl ₃	0.05 g
Seawater.....	500.0 ml
Distilled water.....	500.0 ml

Autoclave at 121C for 15 minutes.