

**ATCC medium: 1951 LB medium (ATCC medium 1065) with 50 mcg/ml ampicillin and 15 mcg/ml tetracycline**

ATCC Medium 1065 (see below) with 50.0 mcg/ml ampicillin and 15.0 mcg/ml tetracycline

*ATCC Medium 1065:*

Tryptone (BD 211705).....	10.0 g
Yeast Extract (BD 212750).....	5.0 g
NaCl.....	10.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.