

ATCC Medium: 200 YM Medium (Agar or Broth)

Agar Medium

YM Agar (BD 271210).....41 g
DI Water.....1000 ml

Autoclave at 121°C.

Broth Medium

YM Broth (BD 271120).....21.0 g
DI Water.....1000 ml

Autoclave at 121°C.

*This product is commercially available from BD.

To make medium from scratch, follow formulation below:

Yeast Extract.....3.0 g
Malt Extract.....3.0 g
Dextrose.....10.0g
Peptone.....5.0 g
Agar (if required).....20.0 g
DI Water.....1000 ml

pH 6.2 +/- 0.2