

**ATCC medium: 2018 *Microlunatus* medium**

Glucose.....0.5 g  
Peptone.....0.5 g  
Monosodium glutamate.....0.5 g  
Yeast extract.....0.5 g  
KH<sub>2</sub>PO<sub>4</sub> .....0.44 g  
(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>.....0.1 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....0.1 g  
Distilled water.....1.0 L

Adjust medium for final pH of 7.0 with NaOH. Autoclave at 121C for 15 minutes.