

**ATCC medium: 209 Marine *Cytophaga* medium**

Nutrient Broth (BD 234000).....8.0 g  
Yeast extract.....5.0 g  
NaCl.....12.86 g  
MgCl<sub>2</sub> . 6H<sub>2</sub>O .....5.29 g  
KCl.....0.75 g  
CaCl<sub>2</sub> . 2H<sub>2</sub>O .....0.74 g  
Fe(NH<sub>4</sub>)<sub>2</sub>(SO<sub>4</sub>)<sub>2</sub> . 6H<sub>2</sub>O .....0.39 g  
Agar (if needed).....15.0 g  
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.